

WATSON®

Dispo Pipette

NEW disposable pipette (5ml/10ml)

Experience high precision graduation achieved by our injection molding technology



Readability of the graduation enhances as aspirating a solution as the contrast increases.



P/No.	Name	Packaging
1650-005S	5ml dispo pipette (Sterilized)	150 pcs.(individual packaging)
1650-010S	10ml dispo pipette (Sterilized)	100 pcs.(individual packaging)

Short scale,
for the ease of
use inside
a clean bench

Body: PS / Adapter:PP
DNase/RNase free
Pyrogen test results:
Endotoxin level less
than 0.05EU/ml
*according to our internal test

<http://watsonbiolab.com>

For Research

Made in JAPAN

Dispo Pipette



Something Different.
WATSON BIO LAB
MADE IN JAPAN SINCE 1988

<http://www.watsonbiolab.com>
info@watsonbiolab.com

Manufactured by:

FUKAEKASEI CO., LTD.

Head Office: 2-2-7 Murotani, Nishi-ku, Kobe
651-2241, JAPAN
TEL +81-78-991-4477
FAX +81-78-991-4491

Sold by:

Watson Bio Lab USA

6540 Lusk Blvd. Suite C158A, San Diego, CA 92121
TEL 949-395-8460

※The specification and the price etc. of the product might change without a previous notice.