


# Pipette Tips



Full size  
 Graduated



**Pipette tip compatibility**

- Watson NEXTY 2-20µl/10-100µl/20-200µl(M)
- Eppendorf Reference 2-20 µl/10-100 µl(M)/50(20)-200 µl
- Gilson Pipetman P20/P100/P200 • Finnpipette Digital 2-20 µl/5-50 µl(M)/10-100 µl/20-200 µl
- Biohit Proline 2-20 µl/5-50 µl(M)/10-100 µl/20-200 µl
- Nichiryo Nichipet EX 2-20 µl/10-100 µl/20-200 µl/5-50 µl(M)/40-200 µl(M) • A&D MPA-200

## Complete & Accurate liquid dispensing! Try our products!

Super thin walled molding for an extra fine edge on the tip.  
 This insures complete dispensing of the liquid in the pipette.



Because the edge is made in the extra fine by the super-thin-walled molding, a more advanced drying well has been achieved.

Cat.No.	Item	Packaging	Unit	Case
110-705Y	200ul, Graduated, Yellow	Bulk	1000tips	10units

### ECO PACK



Just remove and put on. It's very simple!

\*Empty rack Cat#1298-755RA can be used.

Cat.No.	Item	Packaging	Unit	Case
123R-755YS	200ul, Graduated, Yellow, ECO PACK	Stack Rack <b>Sterilized</b>	96tipsX5platesX10stacks	5units

### NEXTY PACK



Thin wall design reduces waste compared to conventional product.

Place the Support Post\* inside the NEXTY Rack(S size) before use. (\*Originally included)

Easy to Refill -just place a refill plate on the NEXTY Rack.

\*For NEXTY RACK use only.

Cat.No.	Item	Packaging	Unit	Case
123NT-755Y	200ul, Graduated, Yellow, NEXTY PACK	Stack Rack	96tipsX5platesX10stacks	5units

### PP RACK




\*This item cannot be used for Nexty Pack and Eco Pack.

Cat.No.	Item	Packaging	Unit	Case
1201-705Y	200ul, Graduated, Yellow	Racked(PP)	96tipsX10racks	5units
1201-705YS	200ul, Graduated, Yellow	Racked(PP) <b>Sterilized</b>	96tipsX10racks	5units



#### System Rack (PC)



The body and lid, made of polycarbonate(PC), does not lose its shape over time use. It can also be autoclaved repetitively.


#### System Rack (PP)



The body and lid, made of Polypropylene(PP) is lightweight, durable, and highly resistant to chemicals and heat. Because of its inexpensive cost compared with System Rack (PC), this System Rack(PP) can be disposed as one time use if used with non-sterilized long tip, or can be reused even after autoclaving if tips are manually refilled.

**Note:** The rack made of PP shrinks slightly after autoclaving, consequently refill plate cannot be used.

#### Refill Plate



One-touch replacement for the System Rack. The paper cover can be recycled.


**Note:** Use a PC (polycarbonate) rack when refilling tips.



Remove the old plate.



Place the new plate. Remove paper skirt.



Refilling tips is complete.





