

Standard Filter Tips




10 μ L
Short


Full size
Graduated



System Rack(PC) System Rack(PP)



Refill Plate



Pipette tip compatibility

- Watson NEXTY/NEXTY-S 0.2-2 μ L/1-10 μ L(Multi)/10 μ L(Fixed volume)
- Eppendorf Reference/ResearchV 0.1-2.5 μ L/0.5-10 μ L(Multi)
- Gilson Pipetman P2/P10 • Finnpiptette Digital 0.2-2.0 μ L/0.5-10 μ L(Multi)
- Biohit Proline 0.5-10 μ L(Multi) • Nichiryo Nichipet EX 0.1-2 μ L/0.5-10 μ L(Multi)
- A&D MPA-10/MPA-20


Cat.No.	Item	Packaging	Unit	Case
1251-204CS	10 μ L, Short Filter Tip, Graduated	System Rack(PC)	Sterilized	96tipsX10racks 5units
1251P-204CS	10 μ L, Short Filter Tip, Graduated	System Rack(PP)	Sterilized	96tipsX10racks 5units
1271-204CS	10 μ L, Short Filter Tip, Graduated	Refill Plate	Sterilized	96tipsX10plates 5units

10 μ L
Long


Full size
Graduated



System Rack(PC) System Rack(PP)



Refill Plate



Pipette tip compatibility

- Watson NEXTY/NEXTY-S 0.2-2 μ L/1-10 μ L(Multi)/10 μ L(Fixed volume)
- Eppendorf Reference/ResearchV 0.5-10 μ L(M)
- Eppen Reference 2-20 μ L(Gray Button) • Gilson Pipetman P-2/P-10
- Finnpiptette Digital 0.2-2.0 μ L/0.5-10 μ L(Multi) • Biohit Proline 0.1-2.5 μ L/0.5-10 μ L(Multi)
- Nichiryo Nichipet EX 0.1-2 μ L/0.5-10 μ L(Multi) • A&D MPA-10/MPA-20


Cat.No.	Item	Packaging	Unit	Case
1252-207CS	10 μ L, Long Filter Tip, Graduated	System Rack(PC)	Sterilized	96tipsX10racks 5units
1252P-207CS	10 μ L, Long Filter Tip, Graduated	System Rack(PP)	Sterilized	96tipsX10racks 5units
1272-207CS	10 μ L, Long Filter Tip, Graduated	Refill Plate	Sterilized	96tipsX10plates 5units

20 μ L


Full size



System Rack(PC) System Rack(PP)



Refill Plate



Pipette tip compatibility

- Watson NEXTY/NEXTY-S 2-20 μ L/20 μ L(Fixed volume)
- Eppendorf Reference/ResearchV 2-20 μ L(Yellow Button)
- Gilson Pipetman P20 • Finnpiptette Digital 2-20 μ L
- Biohit Proline 2-20 μ L • A&D MPA-200


Cat.No.	Item	Packaging	Unit	Case
124-20S	20 μ L, Filter Tip	System Rack(PC)	Sterilized	96tipsX10racks 5units
124P-20S	20 μ L, Filter Tip	System Rack(PP)	Sterilized	96tipsX10racks 5units
126-20S	20 μ L, Filter Tip	Refill Plate	Sterilized	96tipsX10plates 5units

100 μ L

Full size



System Rack(PC) System Rack(PP)



Refill Plate



Pipette tip compatibility

- Watson NEXTY/NEXTY-S 10-100 μ L(Multi)/100 μ L(Fixed volume)
- Eppendorf Reference/ResearchV 10-100 μ L(Multi)
- Gilson Pipetman P-100 • Finnpiptette Digital 5-50 μ L/10-100 μ L
- Biohit Proline 5-50 μ L(Multi)/10-100 μ L • A&D MPA-200

Cat.No.	Item	Packaging	Unit	Case
124-100S	100 μ L, Filter Tip	System Rack(PC)	Sterilized	96tipsX10racks 5units
124P-100S	100 μ L, Filter Tip	System Rack(PP)	Sterilized	96tipsX10racks 5units
126-100S	100 μ L, Filter Tip	Refill Plate	Sterilized	96tipsX10plates 5units



Nuclease-free certificate.

This mark shows RNase/DNase/Human-DNA free. Products with this mark are manufactured in the clean room under our high quality control, and inspected to meet our strict quality standards. RNase/DNase-free certificates can be issued for each lot upon your request.




Full size



System Rack(PC)

Refill Plate



53mm

Pipette tip compatibility

- Watson NEXTY/NEXTY-S 20-200 µL(Multi)/200 µL(Fixed volume)
- Eppendorf Reference/ResearchV 50(20)-200 µL
- Gilson Pipetman P-200
- Biohit Proline 20-200 µL
- Finnpipette Digital 20-200 µL
- Nichiryō Nichipet EX 20-200 µL
- A&D MPA-200

Cat.No.	Item	Packaging	Unit	Case
124-200S	200 µL, Filter Tip	System Rack(PC)	Sterilized	96tipsX10racks
126-200S	200 µL, Filter Tip	Refill Plate	Sterilized	96tipsX10plates






Full size



System Rack(PC)

System Rack(PP)



Refill Plate



53mm

Pipette tip compatibility

- Watson NEXTY/NEXTY-S 20-200 µL(Multi)/200 µL(Fixed volume)
- Eppendorf Reference/ResearchV 50(20)-200 µL
- Gilson Pipetman P-200
- Biohit Proline 20-200 µL
- A&D MPA-200
- Finnpipette Digital 20-200 µL
- Nichiryō Nichipet EX 20-200 µL/40-200 µL(Multi)

Cat.No.	Item	Packaging	Unit	Case
1252-703CS	200 µL, Filter Tip	System Rack(PC)	Sterilized	96tipsX10racks
1252P-703CS	200 µL, Filter Tip	System Rack(PP)	Sterilized	96tipsX10racks
1272-703CS	200 µL, Filter Tip	Refill Plate	Sterilized	96tipsX10plates





Full size

Graduated



System Rack(PC)

System Rack(PP)



Refill Plate



89mm

Pipette tip compatibility

- Watson NEXTY/NEXTY-S 20-200 µL(Multi)/200 µL(Fixed volume)
- Eppendorf Reference/ResearchV 20-200 µL
- Gilson Pipetman P-200
- Finnpipette Digital 20-200 µL
- Biohit Proline 20-200 µL
- Nichiryō Nichipet EX 20-200 µL
- A&D MPA-200

Cat.No.	Item	Packaging	Unit	Case
1252-801CS	200 µL, Long Filter Tip, Graduated	System Rack(PC)	Sterilized	96tipsX10racks
1252P-801CS	200 µL, Long Filter Tip, Graduated	System Rack(PP)	Sterilized	96tipsX10racks
1272-801CS	200 µL, Long Filter Tip, Graduated	Refill Plate	Sterilized	96tipsX10plates






Full size

Graduated




System Rack(PC)

System Rack(PP)



Refill Plate



80mm

Pipette tip compatibility

- Watson NEXTY/NEXTY-S 100-1000 µL/1000 µL(Fixed volume)
- Eppendorf Reference/Research V 100-1000 µL
- Gilson Pipetman P-1000
- Finnpipette Digital 100-1000 µL
- Biohit Proline 100-1000 µL
- Nichiryō Nichipet 100-1000 µL
- A&D MPA-1200

Cat.No.	Item	Packaging	Unit	Case
124-1000S	1000 µL, Long Filter Tip, Graduated	System Rack(PC)	Sterilized	96tipsX10racks
124P-1000S	1000 µL, Long Filter Tip, Graduated	System Rack(PP)	Sterilized	96tipsX10racks
126-1000S	1000 µL, Long Filter Tip, Graduated	Refill Plate	Sterilized	96tipsX10plates





Full size

Graduated



System Rack(PC)

Refill Plate



130mm

Pipette tip compatibility

- Watson NEXTY/NEXTY-S 100-1000 µL/1000 µL(Fixed volume)
- Gilson Pipetman P-1000
- Nichiryō Nichipet EX 100-1000 µL
- A&D MPA-1200
- Eppendorf Reference/Research V 100-1000 µL
- Biohit Proline 100-1000 µL

Cat.No.	Item	Packaging	Unit	Case
124-1200S	1200 µL, Extra-Long Filter Tip, Graduated	System Rack(PC)	Sterilized	96tipsX10racks
126-1200S	1200 µL, Extra-Long Filter Tip, Graduated	Refill Plate	Sterilized	96tipsX10plates





Full size

Graduated



System Rack(PC)

Refill Plate



130mm

Pipette tip compatibility

- Watson NEXTY/NEXTY-S 500-5000 µL/5000 µL(Fixed volume)
- Gilson Pipetman P-5000
- Socorex Acura 0.5-5mL

Cat.No.	Item	Packaging	Unit	Case
124-5000S	5mL, Filter Tip, Graduated	System Rack(PC)	Sterilized	24tipsX1racks
126-5000S	5mL, Filter Tip, Graduated	Refill Plate	Sterilized	24tipsX1plates