

# Device Compatibility Table

\* Please check with the actual product sample if our PCR plates can be used after referring to this acceptable table.

## 137-675C (Half-Skirt)

PCR	
Applied Biosystems	Veriti

Real Time PCR	
Applied Biosystems	StepOnePlus
	7500Fast
	7900HT
	ViiA 7
	QuantStudio 7Flex
	QuantStudio 12KFlex

## 137-674C (Non-Skirt)

PCR	
BioRad	DNA Engine
	C1000-T
	S1000
	MJ Mini (48well)

Real Time PCR	
BioRad	Chromo 4
	DNA Engine Opticon 2
	CFX96 Touch
	Mini Opticon (48well)
Applied Biosystems	StepOne (48well)