

# WATSON®

# Dispo Pipette

## NEW disposable pipette (5ml/10ml)

Experience high precision graduation achieved by our injection molding technology



Readability of the graduation enhances as aspirating a solution as the contrast increases.



Short scale,  
for the ease of  
use inside  
a clean bench

Body: PS / Adapter:PP  
DNase/RNase free  
Pyrogen test results:  
Endotoxin level less  
than 0.05EU/ml  
\*according to our internal test

P/No.	Name	Packaging
1650-005S	5ml dispo pipette (Sterilized)	150 pcs.(individual packaging)
1650-010S	10ml dispo pipette (Sterilized)	100 pcs.(individual packaging)

<https://watsonbiolab.com>

For Research

Made in JAPAN

# Dispo Pipette



Readability of the graduation enhances as aspirating a solution as the contrast increases.



Easy to see aspiration and discharging volume with 2 graduations.



Package design is space saving and easy to stack.

Something Different.  
**WATSON**® BIO LAB  
MADE IN JAPAN SINCE 1988

Manufactured and sold by:

**FUKAEKASEI** CO., LTD.

Head Office: 2-2-7 Murotani, Nishi-ku, Kobe

651-2241, JAPAN

TEL +81-78-991-4477

FAX +81-78-991-4491

E-mail: [info@watson.co.jp](mailto:info@watson.co.jp)

<https://watsonbiolab.com>

※The specification and the price etc. of the product might change without a previous notice.