



Serological Pipettes

Dispo Pipette



Dispo Pipette

High precision scale realized by our injection molding technology.

High precision scales



Precise & easily-recognizable scale realized by our advanced injection molding technology.

2 kinds of scales



Easy to see aspirated or dispensed volume with 2 different kinds of scales employed.



5mL
1650-0055

10mL
1650-0105

Length : 224 mm
Real scale

Length : 244 mm
Real scale

Each adaptor is made of PP and contains PE filter.

Short length

Easy-to-handle in a clean bench with limited handling space.
Length: 224 mm for 5 mL and 244 mm for 10 mL.

Clean

STERILITY ASSURANCE LEVEL (SAL) 10⁻⁶

PYROGEN FREE
Less than 0.05 EU/mL

Sterilized, RNase/DNase-free
Endotoxin level: less than 0.05 EU/mL
Source: Internal inspection.



Space-saving and easy-to-stack package design.



Material: Polystyrene (body) Polypropylene (adapter) Polyethylene (filter)

Tip outer diameter Φ 3.2 mm

Cat. No.	Item	Unit	Case
1650-0055	5 mL Disposable Serological Pipette Sterilized	150 pieces (individual packaging)	5 units
1650-0105	10 mL Disposable Serological Pipette Sterilized	100 pieces (individual packaging)	5 units

